
Fluorescent Dopamine (DnsylD-1™)

Characteristics:

- Dopamine linked through a 5-carbon linker to dansyl.
- Large spectral shift between excitation and emission maxima. Ex:/Em: 333/515nm
- Minimal surrounding autofluorescence at detection wavelengths
- Stable in physiological buffers
- Emits bright blue fluorescence
- Conjugated with small fluorescent tag that maintains the affinity and selectivity of native dopamine
- Blockable by standard inhibitors to dopamine receptors and transporter
- Binds dopamine receptors and DAT at nanomolar concentrations. Use standard inhibitors and methods to distinguish uptake from receptor binding

Specifications:

Weight/Unit: 0.5 mg (1µmole) 1 mg (2 µmole); 5 mg (10µmole)

MW: 499.62 (Confirmed by MALDI)

Appearance: Lyophilized white solid

Solubility: DMSO

Spectral Characteristics: Ex/Em: 333/515 nm

Use at 10X Kd of unlabeled dopamine for respective dopamine receptor isoform.

Purity: 97% (HPLC)

Storage (solid): 4°C, dark, desiccate. DMSO solubilized material: Aliquots at -20°C

Recommended Directions For Solubilization:

(Other solubilization volumes are applicable depending on the mass of the starting material and on the experimenter's final assay concentration).

1. 0.5 mg DnsylD-1

Add 10 µl of DMSO, vortex and solubilize. Centrifuge and collect supernatant. Aliquot and freeze suspended material at -20°C or -80°C. Dilute aliquot in physiological assay buffer to final assay concentration. As a starting point, assay DnsylD-1 at 10-100X the Kd of unconjugated dopamine to the dopamine receptor subtype of interest. Employ a similar method to determine the concentration used for uptake assays.

2. 1 mg DnsylD-1

Add 10 µl of DMSO, vortex and solubilize. Centrifuge and collect supernatant. Aliquot and freeze suspended material at -20°C or -80°C. Dilute aliquot in physiological assay buffer to final assay concentration. As a starting point, assay DnsylD-1 at 10-100X the Kd of unconjugated dopamine to the dopamine receptor subtype of interest. Employ a similar method to determine the concentration used for uptake assays.

3. 5 mg DnsylD-1

Add 50 µl of DMSO, vortex and solubilize. Centrifuge and collect supernatant. Aliquot and freeze at -20°C or -80°C. Dilute aliquot in physiological assay buffer to final assay concentration. As a starting point, assay DnsylD-1 at 10-100X the Kd of unconjugated dopamine to the dopamine receptor subtype of interest. Employ a similar method to determine the concentration used for uptake assays.