

N-kappaB p65 Sandwich ELISA Kit Contents

Materials provided with the kit	96 determinations	
User manual		1
Closure plate membrane		2
Sealed bags		1
Microelisa stripplate	1	4°C
Standard: 250 – 1800 pg/ml (varies per species*)	0.5mlx1 bottle	4°C
Standard diluent	1.5mlx1 bottle	4°C
HRP-Conjugate reagent	6mlx1 bottle	4°C
Sample diluent	6mlx1 bottle	4°C
Chromogen Solution A	6mlx1 bottle	4°C
Chromogen Solution B	6mlx1 bottle	4°C
Stop Solution	6mlx1 bottle	4°C
Wash solution	(20mlx30 fold)x1 bottle	4°C

Storage: 4°C, keep dark

Safety: Avoid contact and inhalation. Stop solution contains 1N sulfuric acid

* For standard concentration, see documentation with kit.