

## Glass Micropipettes – Pre-Pulled Microcapillaries

Part #	Tip ID ( $\mu\text{m}$ )	Base OD* (mm)	Length* (cm)	Price (30 pipettes)	Description and Further Options In Tip Size
<b>MGM 1A-0.3</b>	<b><math>0.3 \pm 0.2</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1A-1</b>	<b><math>1.0 \pm 0.5</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1A-5</b>	<b><math>5.0 \pm 2</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1B-10</b>	<b><math>10 \pm 5</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1C-20</b>	<b><math>20 \pm 10</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1C-30</b>	<b><math>30 \pm 10</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1C-40</b>	<b><math>40 \pm 10</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1C-50</b>	<b><math>50 \pm 10</math></b>	1 or 1.5	4.5 or 5.5	<b>\$225</b>	Indicate specific tip ID, base OD, length, and with or without internal filament on orders.
<b>MGM 1D</b>	<b>&gt;60</b>	1 or 1.5	4.5 or 5.5	<b>\$295</b>	Currently unavailable. 60 - 100 $\mu\text{m}$ tip ID pipettes.
<b>MGM 1E</b>	<b>Customized</b>	On request	4, 5, 8	<b>\$345</b>	Customized base OD (for example 2 mm base OD), wall thickness and/or length.

Packages with 100 pipettes: \$610

Mixtures of pipettes sizes are also available to order