

---

## Protease Inhibitor Cocktail with EDTA (Water Soluble, 100X)

### Part PI-2

#### Directions:

Suspend lyophilized material in 1 ml dH<sub>2</sub>O to make a 100X solution.  
Store unsuspended lyophilized materials at -20°C. Store aliquots in dH<sub>2</sub>O at -20°C.

**Contents** (at 100X concentration: after suspension in 1 ml dH<sub>2</sub>O)

1. AEBSF; 50mM.
2. Aprotinin; 30 uM.
3. E-64; 1 mM.
4. Bestatin; 1mM.
5. Leupeptin; 1mM.
6. EDTA; 5mM.

**Safety:** Harmful. Avoid ingestion, eye and skin contact. Wear protective clothing, gloves and eye protection.