

Sandwich ELISA Kit Contents

Materials provided with the kit	96 determinations	
User manual		1
Closure plate membrane		2
Sealed bags		1
Microelisa stripplate	1	4
Standard: 250 – 1800 pg/ml (varies per species*)	0.5ml×1 bottle	4
Standard diluent	1.5ml×1 bottle	4
HRP-Conjugate reagent	6ml×1 bottle	4
Sample diluent	6ml×1 bottle	4
Chromogen Solution A	6ml×1 bottle	4
Chromogen Solution B	6ml×1 bottle	4
Stop Solution	6ml×1 bottle	4
Wash solution	(20ml×30 fold)×1 bottle	4

Storage: 4°C, keep dark

Safety: Avoid contact and inhalation. Stop solution contains 1N sulfuric acid

* For standard concentration, see documentation with kit.